

## **Enhanced detection of *esr1* mutations in breast cancer: validation of the galeas upcr:esr1 kit on multiple commercially available reference standards.**

**Abstract Submitter:** Michael Parks, United Kingdom\*

Co-Authors: Agata Stodolna, Charlotte Morgan

\*Nonacus Ltd

### **Abstract**

#### **Background**

Mutations in the *ESR1* gene are a significant cause of resistance to hormone therapy in estrogen receptor (ER)-positive breast cancer patients and occur in 10-50% of metastatic endocrine-therapy resistant cancers. Detecting these mutations through liquid biopsy testing can guide treatment towards more effective targeted therapies, as recommended by ESMO and ASCO guidelines. However, current *ESR1* variant testing on cfDNA is mostly conducted using complex, time-consuming, and expensive sequencing technologies.

#### **Objectives**

To validate the Galeas uPCR:ESR1 kit developed by Nonacus using commercially available reference standards provided by SensID, ZeptoMetrix, and SeraCare; and thus confirming the kit's performance in detecting common *ESR1* variants with high analytical sensitivity and specificity.

#### **Methods**

The Galeas uPCR:ESR1 kit uses proprietary qPCR technology developed at Nonacus to detect 11 of the most common *ESR1* mutations, which account for over 90% of variants associated with endocrine therapy resistance. The analytical sensitivity of the kit is greater than 90% at limits of detection ranging from 0.04% to 0.33% variant allele frequency (VAF) across the 11 targets, with an analytical specificity of  $\geq 98\%$ .

#### **Results**

The technology was evaluated using *ESR1* cfDNA reference material provided by SensID, which included the E380Q, D538G, Y537S, Y537N, S463P, Y537C, L536H, and L536R variants. At 15 ng input, all variants were detected at contrived VAFs of 1%, 0.3%, and 0.1%. Additionally, the kit's performance was evaluated on the Zeptomatrix *ESR1* Control Kit, confirming that 10 out of 11 mutations could be detected down to 0.1% VAF, with only the S463P variant not detected at the lowest VAF tested. Further validation with SeraCare ctDNA *ESR1* Mutation Mix AF1 reference material showed positive amplification for all tested VAFs (1%, 0.3%, and 0.1%), with HEX channel variants E380Q and S463P detected down to 0.1% VAF in all replicates.

#### **Conclusions**

The Galeas uPCR:ESR1 qPCR kit presents a simple yet powerful solution for the detection of *ESR1* mutations linked with endocrine therapy resistance in breast cancer. This new technology offers a faster, accurate, and cost-effective alternative to existing sequencing methods, potentially improving treatment outcomes for breast cancer patients. The validation with SensID, ZeptoMetrix, and SeraCare reference standards further supports the kit's robustness and reliability.

#### **Do you have any conflicts of interest?**

Yes, I have a conflict of interest.

I am employed by Nonacus Ltd